Micro Probe System





Index

*NEXTRON

MICRO PROBE SYSTEM



Index



14 Customers



1



Micro Probe System is suitable to analyze the Electrical & Optical properties of the materials. Its advantage is the in-situ measurement of the electrical and optical properties under the various environmental conditions; Vacuum, Temperature, Gas Flow, Humidity, Light Irradiation and X-ray.

NEXTRON R&D Team is Ready to Discuss Complex Requirements

Overview	Specific requirements should be provided	for further assessment.					
Modeling	A dedicated team collaborates directly with clients to design customized chambers that align with specific requirements. A 3D image copy is provided along with a budgetary quote. If necessary, this process is repeated until all requirements are met.						
Production	Custom orders are manufactured using a to ensure precision and quality.	dvanced high-tech machining equipment					
Custom list	 Dimension Number of Probes Optical Components Temperature Range 	 RF Measurement Humidity Control High Voltage Probe Tip 					
	 Gas (Mass Flow Controller) 	 Additional specifications as 					

 Additional specifications as required



1



RT ~ +1000 °C (Model: MPS - CHU) Sample stage size: Ø 1/2"

RT ~ +750 °C (Model: MPS - CHH) Sample stage size: Ø 1/2"

RT ~ +450 °C (Model: MPS - CHL) Sample stage size: Ø 1/2"

-40 ~ +200 °C (Model: MPS - PT) Sample stage size: 19 x 19 mm

-40 ~ +170 °C (Model: MPS - PTH) Sample stage size: 16 x 16 mm

80 K ~ RT (Model: MPS – LN) Sample stage size: Ø 15mm

Main Probe System Structure

MPS - PTH





1

Quick & Easy

General





Motorized Probe Module

Precise & Positioning under vacuum



Specification

Module Speed: 5 mm /s Stroke (X , Y axis): 8 mm Stroke (Z - axis): 3 mm Holding Force: 150 gf – 200 gf Manipulate Dimension: 20 × 20 × 30 mm Resolution: < 1 µm

The manual probe module is compact, measuring $30 \times 20 \times 20$ mm. It enables rapid and stable contact with contact pads exceeding 500 µm using a manual contact method. Due to its flexible motion, the module smoothly adapts to the displacement of the contact point resulting from thermal expansion, maintaining stable contact.

video at: https://www.youtube.com/watch?v=STMf0xmsdQo&t=2s video at: https://www.youtube.com/watch?v=dx8p2lL6QP4 Ensuring stability and high resolution is crucial in precision applications. The Piezo-Driven Positioner offers an optimal solution by enabling remote micro-positioning. Nextron piezodriven micro probe system allows for precise placement at the desired location. Even when contact pads are extremely small, the system enables precise measurements using a 10 µm probe tip. Additionally, the piezo-driven micro probe system is competitively priced compared to other market alternatives.

video at: https://bit.ly/3Ne5tkQ



Probe Tips

Tungsten Probe



Tungsten Probe Infor	mation
----------------------	--------

Tip Dia.	Plate	Operating temp./time @1atm		
50 µm	None 450 °C / 72 h			
100 µm	None	450 °C / 72 h		
	Special allov	450 °C / 72 h		
	Special alloy	750 °C / 3.5 h		

* Small tip size

Nextron 1K Probe



Nextron 1K Probe Information



* The unique probe capable of withstanding extreme he up to 1000°C in air.

- Tungsten probes with tip diameters of 50 μm and 100 μm are recommended for applications requiring small contact pads. For enhanced precision, the use of an optional precise probing arm is advised.
- · Larger tip sizes provide increased sample clamping force and durability.
- · The table above presents recommended models for appropriate probe selection.
- All probes feature a spring structure within or outside the probe module, ensuring consistent contact stability in gas flow and scanning experiments. A key advantage is that the contact remains stable even in the presence of movement caused by thermal expansion and vibration.

Nextron 1K Probe High-temperature resistance test



The graph illustrates the resistance stability of the Nextron 1k Probe under conditions of 1000 °C at 1 atm air over a 24-hour period.

Equipment Used:

- · MPS-CHU probe system with Nextron 1K Probe
- SMU: Keithley 2450



6 Versions Depending on Temperature

Peltier Stage

2







MPS - PTH

Temp.: -40 ~ +170 °C Sample stage size: 16 × 16 mm Stage material: Special alloy coated Cu Ø 3.2 mm bottom hole



Temp.: MPS - CHL: RT ~ +450 °C MPS - CHH: RT ~ +750 °C

MPS - CHU: RT ~ +1000 °C Sample stage size: Ø 1/2" Stage material: Alumina

6 Versions Depending on Temperature

Liquid Nitrogen Cooling Sample Stage



video at: https://bit.ly/3FjcGv8

* NEXTRON MICRO PROBE SYSTEM

Manual Probe Chambers



6-Channel Micro Probe System



Micro Probe System

8-Channel Micro Probe System



10-Channel Micro Probe System





2

The In-situ Raman Battery Test Cell is considered a breakthrough product in the rechargeable battery industry.



In-situ Raman battery cell is capable of measuring the Raman spectrum of the cross-section of a battery sample. The In-situ Raman battery cell is developed for the in-situ study of the interface between layers such as anode, cathode, and separator. It is available to measure in-situ Raman not only liquid electrolyte batteries but also solid-state batteries. The product has a feature that can easily inject a liquid electrolyte and simple preparation of test samples.

Specification

- Cell dimension: 55(l) × 55(w) × 30(h) mm
- Inner electrolyte dish volume: 150 µL
- Cell material: PEEK
- Corrosion resistance:
- All metal parts are coated with platinum
- Battery charge/discharge holder type: M3 bolt clamping type
- Battery holder material: Stainless 316L

- Weight of cell: 150 g
- Number of electrode: 2 ea
- Cathode/anode electrode type: Banana plug
- Sample total thickness: 300 µm
- Sample loading size: 15 × 6 mm
- Window material: Sapphire
- Window dimension: Ø 25 × t 1.0
- Minimum optional working distance: 1 mm
- Electrolyte inlet/outlet fitting: 1/16 inch Fingertight Fitting













*NEXTRON

MICRO PROBE SYSTEM

3



3

Sample loading to the electrode holder

- Align the sample
- Align the sample on Holder A as required. 1. Soak the separator in the electrolyte prior to use where applicable.





- When using a dish place the dish and the rod into the designated grooves. 2. Ensure the end of the separator is fully contained within the dish.



- * The total thickness of the sample is approximately 300 µm.
- * For samples with significant thickness differences, a holder suitable for the sample size is to be used.

Cover electrode B and tighten the bolt. 3. Raising the dish allows for easier tightening of the bolt.



Cut the protruding sample to fit the cross-sectional side of the electrode holder. 4.









Hall Measurement Micro Probe System



The Hall Measurement Micro Probe System is designed for applications requiring integration with conventional electromagnets, facilitating electric transport experiments, especially Hall effect studies. This system is compatible with Nextron's environmental controllers, enabling comprehensive environmental control including temperature, vacuum, humidity, and gas management. Equipped with Nextron's unique four-probe configuration, this system allows for advanced magnetic in-situ experiments, providing precise measurements and analyses in controlled conditions.

Specification

Min. air gap	40 mm (both in-plane & out-of-plane)
Max. pole face	3 inch
	HMMPS-PT -40 ~ 200°C
Temperature ranges	HMMPS-CH RT ~ 1000°C
	HMMPS-LN 80K ~ RT (coming soon)
Chamber cooling option	Active TEC air cooling block (Required for high temperature model; HMMPS-CH)
Motorized chamber rotation option	Angular magnetoresistance measurements between in-plane & out-of-plane (coming soon)
Chamber mounting stand option	Easily position the chamber between the magnetic poles (Custom made to fit your magnet)
SMIL switching hox option	Useful for measuring Hall voltage
Sivio switching box option	in a van der pauw sample (coming soon)
Number of probes	4 probes
	Vacuum tight < 1 × 10 ⁻⁴ mbar
Environmental control	Humidity control option 4 ~ 95 %RH
Environmental control	Vacuum partial pressure control option 0.001~1300 mbar
	MFC Gas control option
Dimension	150 × 80 × 80 mm
Signal port	SMA (SMA to BNC cable included)
Vacuum & gas port	1/4 inch & 1/8 inch Swagelok Tube fitting
Optical window option	Customizable for optical access
	Leakage Current 100 fA
Electrical characteristics	Maximum DC Voltage/Current 300 VDC/1000 mA
	Frequency DC ~ 300 MHz
	HMMPS-CH Ø 15 mm
Thermal stage size	HMMPS-PT 19 × 19 mm
	HMMPS-LN Ø 14 mm



Micro Probe System for THz Signal **Transmittance Measurement**



for test and reference samples

This probe system is perfectly suited for terahertz transmission experiments.

It features a temperature-controlled sample holder with dual transmission apertures. It allows for simultaneous comparative analysis of test and reference samples.

The unique Nextron probe ensures stable electrical probing and sample fixation, even in vertical orientations.





Specification

- Dimension: 140(l) × 70(w) × 31.6(h) mm
- Stage size: 32 × 18.2 mm
- Stage transmission hole: Ø 2 mm
- Temperature Range:
- 80 K~RT(LN model)/ -30 ~ 150 °C(PT model)
- Temperature Accuracy: ±0.1 °C
- Temp Resolution: ±0.1 °C
- Max. Ramping Rate: 30 °C/min
- · Top and bottom window:
- Ø 43 mm, Fused silica 2T





SPM(scanning probe microscopy) Probe System



Combine precise temperature control with electrical probing for **in-situ electrical and topographical measurements**.

Specification

- Dimension: 90(l) × 35(w) × 20(h) mm
- Stage size: 16 × 16 mm
- Temperature range: -30 ~ 150 °C¹⁾
- Max. Ramping Rate: 30 °C/min
- Temperature Accuracy: ±0.1 °C

¹⁾The lower limit of the temperature range is typically attained when a coolant is approximately 0 °C.

SPM(scanning probe microscopy) Thermal Stage



Specification

- Dimension: 20(l) × 26.5(w) × 6.5(h) mm
- Stage size: 16 × 16 mm
- Temperature range: RT ~ 150 °C
- Max. Ramping Rate: 30 °C/min
- Temperature Accuracy: ±0.1 °C



XRD Micro Probe System (Environment System with Probes for Mounting XRD Machine) 21.5 mm 23.5 mm 0.**.0 Bottom to Stage XRDMPS-CHU 24 mm 75 mm

Properties		Specification	Properties		Specification
Temperature Range		RT ~1000 °C		20	0 ~ 180°
Vacuum	Rotary Pump	~ 10 ⁻³ mbar at RT	Angle	Ψ	0 ~ 85°
Vacuum	Turbo Pump	~ 10 ⁻⁵ mbar at RT		Φ	0 ~ 360°
Heating Stage Size		Φ 1/2"	Dome Material		PEEK
Number	of Probes	4 ea	Weight		450 g

D8 Bruker XRD set up image







Technical Description

The XRD Micro Probe System is a probe station designed for in-situ XRD measurements under high temperature and high vacuum conditions. The system incorporates Nextron's unique manual probe, enabling the application of voltage and the measurement of current. This system is particularly effective for electric-field-induced phase transitions in ferroelectric materials, resistance analysis during PRAM phase changes, and ion-based battery studies.

More Information: Applications-XRD: In-situ Observation of Structural- and Resistance-Changes by Redox Reactions

Applications

The XRD Micro Probe System is applicable to a wide range of research fields, particularly for in-situ analysis within XRD chambers. The system supports the following applications:

- Phase transition analysis under oxidizing or reducing gas atmospheres, utilizing XRD for structural observations and probes for electrical property measurements.
- Crystal structure analysis of materials subjected to periodic temperature variations, such as battery cells.
- In-situ characterization of anode and cathode crystal structures in battery cells during charging and discharging cycles.
- Phase transition analysis during Chemical Vapor Deposition (CVD) processes, particularly for monitoring thin-film thickness variations.
- · Electric-field-assisted phase transition studies under varying temperature conditions.
- Temperature-dependent phase transition analysis, incorporating both XRD measurements and electrical property evaluations.

3



Probe Head Chamber - Pogo pin

3



The probe head can be customized according to sample.





Micro Probe System for LCC Chip



Temperature Range	28-pin	48-pin	64-pin
-30 ~ 150 °C	LCC MPS - PT28	LCC MPS - PT48	LCC MPS - PT64
80K ~ RT	LCC MPS - LN28	LCC MPS - LN48	LCC MPS - LN64

* Customization is avaible on request.





KICRO PROBE SYSTEM

Feature Comparison of Probe Head System and LCC System

Feature Comparison of Probe Head System and LCC System						
Product	Probe Head System	LCC System				
Model Name	MPW-WXXXX	LCCMPS-XX				
Discription	Using Pogo pin	Using LCC socket and Chip Carrier				
Contact Pad size	0.5mm or more (Direct Contact)	Less than 0.5 mm (Wire Bonding)				
Pitch	> 2.5 mm	-				
Remark	Customize based on customer sample's electrode information	28, 48, 64 Pins available				
Temperature	-30 to 150 °C	80K to RT -30 to 150 °C				

MEMO





Micro Probe System for Field Emission Measurements

4 probes FEMPS-PT6C 6 SMA feedthough for probes and electrodes NW16 Slidable top anode Special alloy coated stainless plate Option: Piezo-driven XZ manipulator Cathode sample stage Special alloy coated stainless plate

Optical Window Stage Chamber for Integrated Electro-Optical Experiments



Bottom side





RF Micro Probe System with GSG Probe



GSG Probe Tip



Technical Specification

Standard Pitch Range (µm)	Minimum Pad Size (µm)	Maximum Frequency (GHz)
50 ~ 1250	80 × 80	40

Temperature Range

RFMPS-VAC2C	RFMPS-LN2C	RFMPS-PTH2C	RFMPS-PT2C	RFMPS-CHL2C
Non-temp. control	80K ~ RT	-40 ~ 170°C	-40 ~ 200°C	RT ~ 450°C



G

Micro Probe System with Wire Push Contact

MPS-PT6C

Micro Probe System without Thermal Stage





Micro Probe System with Terminal Block

3

6 Probes and 24-Pin Terminal Block



43

Micro Probe System for Thermoelectric Devices (Dual Peltier Model)









Micro Probe System

13-channel Micro Probe System



Ø 162 mm





Micro Optical Chamber (MOC)

MOC is miniature optical chamber (MOC) for optical measurements.



Low weight (300 g), small volume (< 30 cc), low profile (24.5 mm)
Focal distance from sample stage to top of window: ~3 mm
6 versions depending on various temperature ranges

-40 ~ +200 °C (Peltier type, MOC-PT)
 -40 ~ +170 °C (Peltier type, MOC-PTH)
 RT ~ +450 °C (Ceramic Heater type, MOC-CHL)
 RT ~ +750 °C (Ceramic Heater type, MOC-CHH)
 RT ~ +1000 °C (Ceramic Heater type, MOC-CHU)
 80 K ~ RT (Liquid Nitrogen type, MOC-LN)

0...0

 \bigcirc



Piezo Module Micro Probe System

KICRO PROBE SYSTEM



For the contact pad smaller than 400µm x 400µm, Piezo type probe is recommended.

Program		Ye Probe A.
N-PZM	– 🗆 X	
Probe A	Probe C	
y Avia	v Avia	
Voltage 25.0 EV1	Voltage 25.0 EVI	
Frequency (60.0 Diffs)	Frequency 50.0 [kHz]	
Duration 100,0 (puter)	Duration 100,0 [pulse]	-
Detault Save	Default Save	0.4
y Axia	y Axis	
Voltage 25.0 [V]	Voltage 25.0 [V]	
Frequency (60,0 (kHz)	Frequency: 50,0 [kHz]	
Duration 100,0 [puise]	Duration: 100,0 [puise]	
Default Save	Default Save	
z Axia	z Axiz	
Voltage [2]	Voltage 25,0 [V]	
Frequency (60,0 [kHz]	Frequency: 60,0 [kHz]	
Duration 100.0 [puise]	Duration 100.0 [puise]	
Default Save	Default Save	
and a D	Dank a D	
TUDE D	PTODE D	
E AXIE Voltage (TE D	X AXIX	
Formage (210 Phil)	Francisco FRA	
Duration (100.0 Found)	Duration 1010 [prize]	
Default Save	Detault Save	
v Aziz	v Ariz	
Voltage (20,0 [V]	Voltage 25,0 [V]	
Frequency (60,0 [kHz]	Frequency 60,0 [kHz]	
Duration 100.0 [puise]	Duration 100.0 [pulse]	
Detault Save	Default Save	1
z Axis	z Axis	
Voltage (25,0 [V]	Vollage 25,0 [V]	
Frequency 60.0 BHz]	Frequency: 60.0 [kH2]	
Received and Received	Duration: 1010 [pulse]	
Doueupou innnn ibnieer	Portanett lowe aportees	



- D X Probe A



4

Dimension



(Unit: mm)

			w	h ¹⁾	a ²⁾		f ³⁾	р	
Micro Probe System	4-ch	140	70	30.3	-	11.3	27	5.3	6
	6-ch	140	80	30.3	-	11.3	27	5.3	6
	8-ch	140	99	30.3	-	11.3	27	5.3	6
	10-ch	180	99	30.3	-	11.3	27	5.3	6
	13-ch	-	-	30.3	160	11.3	27	5.3	6
Piezo Micro Probe System	4-ch	180	99	33.3	-	11.3	27	5.3	6
	8-ch	-	-	33.3	220	11.3	27	5.3	6

¹⁾ With the addition of a cooling plate, the CHH/CHU model's height increases by 6.5mm.

²⁾ Apothem in a polygon shaped chamber (Distance between Opposite Sides)

³⁾ Standard 1/4" swagelok tube fitting, 23.2 mm with 1/8" tube fitting, 33 mm with NW16 option.

Thread Information for Mounting and Bracket Design



(Unit: mm)

		I	w
Micro	4-ch,6-ch	125	50
Probe System	8-ch	125	75
Piezo Micro Probe System	1 ab	168	85
XRD	4-cn	64	38
RBC	46	46	

In-situ Stretchable Probe System



Our system Provides in-situ electrical characterization of planar and flexible devices during stretching, eliminating the limitations of conventional clamp-based methods.

video at: https://youtu.be/Loh6-39Y5q0?si=H4CYmf6ODmcl477W

Specification

- Effective stroke: 10 ~ 30 mm
- Max sample width: 30 mm
- Max stretching speed: 10 mm/s
- Cycle repeatability: 10 μm

- Max tensile force: 3 kg
- Force meter resolution: 5 gf
- Electrical measurements:
 2 or 4 wire measurements (500 VDC/1 A)
- Probing contact force: 5 ~ 20 gf

The Stretchable Probe System fits a periodic expansion/return experiment with an electrode coated on a stretchable substrate.

The probe, placed directly on the electrode located on the stretchable substrate between the clamps, allows for studying changes in electrical characteristics in various environments such as vacuum, gas, humidity, and UV.

In a typical experiment, electrical resistance is measured by securing the ends of a flexible substrate with clamps and connecting wires to these clamps. However, the clamping force exerted on the sample during stretching can introduce variability, making it difficult to obtain reliable data.

With the Stretchable Probe System, by positioning the probe directly onto the center of the sample, the variable introduced by clamp pressure, a common issue in traditional methods, is completely eliminated. Consequently, this approach facilitates the analysis of intrinsic properties related to adhesion between the elastic substrate and the electrode.

Furthermore, it is expected that the characteristics of FET devices implemented on the stretchable substrate can be measured in real-time, in accordance with the expansion/ return of the stretchable substrate.





Gas Sensor Measurement System (GSS)



Optimized system for analysis according to temperature and gas atmosphere. Especially gas sensors.

Specification

- Small size for compact set-up: 324(l) × 233(w) × 170(h) mm
- Two options of control range: GSS-PT(10 ~ 200 °C) / GSS-CHL(RT ~ 450 °C)
- High accuracy of mass flow control (±0.4 % of reading and ±0.2 % of full scale)
- Various gas selection: 150+ gases including corrosive ones & 20 settings of mixed gas

Gas Selection

Mass flow	Single	gas ¹⁾	Mixed gas
101035 11000	Standard	Corrosive	20 definable das mixes
500 SCCM	128 selectable gas	32 selectable gas	20 definable gas mixes

¹⁾ Get more information about gas selection on the contact page of "microprobesystem.com"

Display for Temp. & Gas flow



Test chamber

Select between PT and CH models according to temperature range



Highly Configurable and Adaptable System



Integrated Software

Software Main Display



Simply set your recipe and execute the process!

Nextron's integrated software enables seamless control of Keithley 2400/2450 as well as various environmental control systems. The software facilitates convenient and efficient management of temperature, humidity, gas flow, pressure, and electrical measurements within a single platform.

NEXTRON



System Configuration for Gas Sensor Analysis



MEMO





Vacuum Partial Pressure Control System

Micro Probe System with Baratron Gauge



Vacuum Partial Pressure Controller



Pressure Control Range

0.001 ~ 1000 Torr

The Precision Vacuum Pressure control system can control partial pressure with the desire gas. The software can automatically acquest and maintain it as a recipe.

video at: https://bit.ly/3yiOW91

Precision Humidity Control System





Humidity sensor¹⁾ can be optionally installed inside the all chambers of Micro Probe System.

The Precision Humidity Control System consists of a two-channel MFC. DI water

Precise control of relative humidity can be

quickly achieved PID control method.

bubbling bath and humidity control

software.



Technical Specification

Typical RH Accuracy (%RH)	RH Resolution (%RH)	Normal Ramp Speed (%RH / min.)	Nomal Range
2 (@1 ~ 97 %RH)	0.01	10	4 ~ 95 %RH ²⁾

¹⁾ Fully calibrated, linearized, and temperature compensated output.

²⁾ Ranges may differ depending on experimentation and laboratory temperature.



8

Controllers

Precision MFC Gas Control Stage

Gas Flow Control



The Ultimate Solution for Precise and Stable Flow Control

The PMGC -XX Gas Flow Station, built with the mass flow control unit, gas flow pipeline, and electrical components, is the ideal solution for various applications demanding precise and stable flow control.

Versatile Gas Control

The Precision MFC gas control stage offers exceptional flexibility, supporting over 130 userselectable gases. Whereas conventional MFCs are restricted to a single gas type, this system enables seamless gas switching and precise concentration adjustments. Its compact design ensures efficient integration into various applications requiring dynamic gas control.

Ease of Use

Ready for operation with just a few gas lines and a USB connection.

Comparison of MFC Gas Control Stage

Nextron MFC Stage

Typical Competitors





	NEXTRON	Typical Competitors	Remarks
Steady State Control	0.01% ~100%	5% ~100% (<mark>20:1</mark> turndown ratio)	Low full-scale
Range	ratio)*	2% ~ 100% (<mark>50:1</mark> turndown ratio)	Large full-scale (>500sccm)
Flexibility of Gas Connections	130+ user-selectable gases	Only a single gas type with a fixed concentration	Refer to diagram

*It has a turndown ratio of 100:1, even for corrosive gases, which is higher than that of typical competitors.



TEM Sample Cleaning / Gas Reaction Chamber



Technical Specification

Maximum Temperature	Max. Heating / Cooling Rate
1000 °C	1 °C/sec

Stage for Eclipse Ti-U Inverted Microscopes



Upper Side

Under Side







Polymer Capacitance Measurement Fixture

Polymer Fixture Equipped on the Peltier Stage





The polymer fixture is designed for measuring the dielectric properties of polymer films. Using this fixture, characteristics such as temperature-capacitance and frequency-capacitance can be analyzed in various environments, including vacuum, humidity, and gas conditions.

The reliability of measurements is enhanced as the fixture ensures consistent contact pressure using a single probe module. Additionally, a guard ring generates a uniform electric field, improving measurement precision. The fixture also provides well-defined top electrodes, facilitating accurate and reproducible measurements.

 This fixture is not an off-the-shelf product and is supplied through customized manufacturing after consultation with customers.

video at: https://youtu.be/N31PeDQ7AEM







Options

Optical Covers

ZnSe Window for IR Transmittance

WC-39





Optical Covers

Cover for Microscope Objective Turret



9

*NEXTRON

MICRO PROBE SYSTEM



Options

Optical Covers

Optical Cover Information

		Window S	Selection Guid	e		
Clear Aperture (mm)		18			43	
Thickness ¹⁾ (mm)	1.0±0.1	0.50±0.05	1.0±0.1	2.0±0.1	1.00±0.05	2.0±0.1
Material ²⁾	Fused Silica	Sapphire	Zinc Selenide	Fused Silica	Sapphire	Zinc Selenide
Wavelength ³⁾ (nm)	200 - 2200	330 - 5500	600 - 18000	200 - 2200	330 - 5500	600 - 18000
TWD ⁴⁾	λ/10	2λ(typical)	λ/10 @10.6µm	λ/10	2λ(typical)	λ/10 @10.6µm
Surface Quality ⁵⁾	20-10	80-50	60-20	20-10	80-50	60-20
Structure	amorphous	Random Axis Orientation ⁶⁾	Isotropic (Cubic Zincblende)	amorphous	Random Axis Orientation ⁶⁾	Isotropic (Cubic Zincblende)
WD (mm)	7±0.2	6.5±0.15	7±0.2	8±0.2	7±0.15	8±0.2

* In addition to the standard options listed in the table, windows can be customized to meet customer specifications, including size, coating, material, etc,.

¹⁾Sapphire windows, owing to their higher modulus of rupture, allow thinner vacuum window applications compared to fused silica or ZnSe. The listed thicknesses are the minimum recommended for vacuum use.

²⁾ Micro probe system(MPS) includes a fused silica window as a standard feature at no additional cost.
³⁾ Anti-reflection coating is available for UV, VIS, and NIR wavelengths.

⁴⁾ For sapphire windows, this value(Transmitted Wavefront) represents the Surface Flatness(P-V).

⁵⁾ Scratch-Dig numbers for surface quality conform to the US government standards MIL-PRF-13830B.

⁶⁾ A wafer-grade sapphire window with a C-plane(0001) surface orientation can be selected as an option.

Optical Covers

Windows Transmission Data for Optical Measurements





Options

Signal Cable

Standard: Coaxial Cable(SMA-BNC)







• Length: 1.5 m

Length: 1.5 m
Impedance: 50 Ω
Leakage current: 100 pA

- Impedance: 50 Ω
- Leakage current: 100 fA

Breadboard & Digital Microscope (Optional)



Model	DM-5MP
Resolution	5M pixels (MJPEG codec 2592 × 1944 resolution)
Magnification	10 ~ 140X
Sensor	Color CMOS
Frame rate	10fps @5MP/3MP/2MP 25fps @130M pixels 30fps @VGA
Lighting	8 LED (on/off)
Measuring function	Available
Calibration function	Available
Weight	125 G
Dimensions	10.5(l) × 3.2(w) cm



Electric Characteristics

- Leakage Current: 100 pA standard SMA-BNC coaxial cable
- Leakage Current: 100 fA option SMA-TNC triaxial cable tested by using Keithley 4200-SCS Semiconductor Characterization System.
- Maximum DC Voltage / Current: 500 VDC / 1000 mA



- Frequency: DC ~ 300 MHz
- Measured by using Agilent E5071C network analyzer



Temperature Resolution and Accuracy

Model: MPS-PT



- Temperature Setting Conditions: 0.1 $^\circ \rm C$ heating condition from 30.04 $^\circ \rm C$ to 30.14 $^\circ \rm C$ for 1 hour
- Featuring Peltier thermoelectric elements, the MPS-PT model offers unparalleled temperature stability and control.
- The accompanying graph demonstrates the system's capability to control the sample stage temperature with a precision of ± 0.01 °C.

*NEXTRON

MICRO PROBE SYSTEM



10

Temperature Cycle Test

Temperature cycle test can be performed accurately with a Peltier type Micro Probe System (model: MPS-PT or MPS-PTH).

Process Value

Set Value

video at: <u>https://youtu.be/cUdyocmzVM8</u>

Software





Model: MPS-PT





10

*NEXTRON

MICRO PROBE SYSTEM

Technical Data

Vacuum Test

Model: MPS-CH



Vacuum Leak Test

Model: MPS-PT







Technical Data

Humidity Test

Model: MPS-CH



This Figure demonstrates the precise control of humidity levels between 10 % RH and 90 % RH.

- Rapid response time minimizes settling time and ensures quick adjustments to maintain desired humidity levels.
- Adjustable total flow rate from 100sccm to 1000sccm allows for fine-tuning to match specific application requirements.
- Integration with Nextron Gas Control System enables simultaneous control of humidity and gas concentration.

Heating Stage Temperature (Internal & Surface)

Model: MPS-CHH





10



10



Thermal Isolation

Model: MPS-CH

The MPS chamber is designed with optimized insulation, ensuring that the chamber temperature remains at a low level even during high-temperature operation.



• The measured temperature at the bottom-center of the chamber, the closest spot to the heater, confirms the effectiveness of thermal insulation.

Window Temperature during High-temperature Experiment



Careful attention is required to protect the objective lens of a microscope during hightemperature experiments. The graph provides data to assist in selecting the material for the optical window of the cover.

You can check the temperature of the window according to the absorption of radiant heat by material and the vacuum level in the chamber.



10



Technical Data

X-ray transmittance of PEEK dome

Model: XRDMPS-CHU



Results

- Transmittance (PEEK dome 1.0 t) : ~15 %
- Transmittance (PEEK dome 0.7 t) : ~45 %

Test conditions

- XRD : Bruker D8 Discover
- Power : 40 kV, 40 mA
- Source : Cu Ka1, 1.5046 Å

MEMO





XRD: In-situ Observation of Structural- and Resistance-Changes by Redox Reactions



- XRD chamber: Nextron XRD Micro Probe System
- Sample: SrFe_{0.8}Co_{0.2}O_{2.5} on 001 SrTiO₃
- Beam line: Pohang Accelerator Laboratory 3D beam line
- Measurement: Real-time XRD and resistance measurements in supply of $N_2(2.5~hrs),\,O_2(6~hrs),\,and\,N_2$ (2.5 hrs) in sequence. Flow rate was 20 sccm.
- Results: The lattice constant of SrFe_{0.8}Co_{0.2}O_{2.5} thin film shifted sequentially to 3.99 Å(N₂), 3.90 Å(O₂), and 4.00 Å(N₂). Associated resistance changes is observed due to topotactic transformations.

This research work is supported by 'Busan Open Laboratory Business Meeting Market Demands' project.

Applications



Materials and Design 182 (2019) 107970



Figure. Transmittance hysteresis curves at wavelength of 2000 nm for (a) thermal and (b) IPL sintered VO₂ films. The temperature corresponding to themaximumof first-order derivative curves for (c) thermal and (d) IPL sintered VO₂ films.

450

600

MICRO PROBE SYSTEM

86

Hydrogen Gas Sensor ACS Sens. 2021, 6, 4145-4155





Applications

Hydrogen Gas Sensor

ACS Sens. 2021, 6, 4145-4155



Figure. Integrated in-plane microheater's (a) temperature variations in terms of the applied DC bias voltages, (b) repeatability at 180 °C, (c) sensitivity and linearity, and (d) surface temperature analysis by a thermal imaging camera at 8 V DC bias.





Raman Mapping



2D Molybdenum disulfide(MoS₂)

Applications

Raman Spectrometry

Vanadium dioxide(VO₂) nanowire Study of phase transition by temperature





Photocurrent Mapping



KICRO PROBE SYSTEM



Humidity Sensor

Adv. Mater. Technol. 2022, 7, 2100751



Figure. (a) Schematic diagram of the custom-made early detection of pipeline leak etection system consisting of PtSe₂ nanograin-based humidity sensor, (b) precise detection of input humidity pulses by the PtSe₂-based humidity sensor, (c) response of best-in-class reference humidity sensor for comparison.

Applications

Humidity Sensor

Adv. Mater. Technol. 2022, 7, 2100751



Figure. IV characteristics of ϵ -Ga₂O₃ layers when exposed to different concentrations of (a) H₂, (b) NO₂, (c) O₂, and (d) CO.



*NEXTRON

MICRO PROBE SYSTEM

Light-emitting diodes

Nanoscale, 2019, 11, 18444–18448



Figure. Temperature-dependent I–V characteristics for (a) the reference LED and (b) the self-protective LED under pulse operation.

Applications

Piezoresistive Pressure Sensor

Sensors and Actuators A 314 (2020) 112217



The bridge voltage as function of pressure after an offset voltage correction at 100 kPa.



The Wheatstone bridges connections to the SMU; schematically (left) and through microscope (right).



Electromigration test

J. Phys.: Condens. Matter 34 (2022) 175401



Figure (a) Top view of the sample used in electromigration test. (b) A schematic overview of the electromigration measurement setup.

Applications

Electromigration test

J. Phys.: Condens. Matter 34 (2022) 175401



Figure. Electromigration results for 300 °C-deposited Al stripes under 1 MA cm-2 at 200 °C, 250 °C, 300 °C, and 350 °C at 10 h.





MIT(Metal-Insulator Transition)

Single crystalline VO₂ nanobeams Study of MIT(Metal-Insulator Transition)



Resistive Gas Sensor Characterization



Measured while purging hydrogen and nitrogen mixed gas without using vacuum pump Sample: Palladium coated Zinc oxide($ZnO/A_{12}O_3$) nanowire for hydrogen gas sensing

Applications

Rapid Microfluidic PCR(Polymerase Chain Reaction)

Slide glass substrate Denaturation at 94 °C and annealing/extension at 59 °C



Thermal Resonant Frequency Measurement





11



- Gas sensor characterization
- Rapid Microfluidic PCR(Polymerase Chain Reaction)
- · Photovoltaic cells' and thermoelectric materials' characterization
- Transistors, diodes, LED, ... testing
- Bulk and thin film materials' thermal conductivity measurements(3ω-Method)
- Phase transition materials' electrical/optical characterization (metal oxides, Memristor,...)
- Characterization of MEMS/NEMS mechanical and electro- mechanical resonators(reference clocks, mass sensors)
- Characterization of micro-coils and micro-antenna for inductive sensors(Impedance spectroscopy of biological tissues, In vivo RMN)
- Capacitive, resistive & resonant micro/nano sensors testing





Zhen Cui, Yaqian Zhang, Dong Hu, Sten Vollebregt, Jiajie Fan, Xuejun Fan and Guoqi Zhang "Effects of temperature and grain size on diffusivity of aluminium: electromigration experiment and molecular dynamic simulation." J. Phys. Condens. Matter 34,175401(2022).

Sihyeok Kim, Gurpreet Singh, Mintaek oh, and Keekeun Lee "An Analysis of a Highly Sensitive and Selective Hydrogen Gas Sensor Based on a 3D Cu-Doped SnO₂ Sensing Material by Efficient Electronic Sensor Interface." *ACS Sens.* 6,4145-4155 (2021).

Paul Grey, Manuel Chapa, Miguel Alexandre, Tiago Mateus, Elvira Fortunato, Rodrigo Martins, Manuel J. Mendes, Luís Pereira "Combining Soft with Hard Condensed Matter for Circular Polarized Light Sensing and Logic Operations." *Adv. Opt. Mater.* 9,2001731(2021).

Nor Syahira Mohd Tombel, Siti Aishah Mohamad Badaruddin, Firzalaila Syarina Md Yakin, Hasan Firdaus Mohd Zaki, Mohd Ismahadi Syono, "Detection of low PPM of volatile organic compounds using nanomaterial functionalized reduced graphene oxide sensor." *AIP Conf Proc* 2368, 20004 (2021).

Taejung Kim, Seungwook Lee, Wootaek Cho, Yeong Min Kwon, Jeong Min Baik and Heungjoo Shin, "Development of a Novel Gas-Sensing Platform Based on a Network of Metal Oxide Nanowire Junctions Formed on a Suspended Carbon Nanomesh Backbone." *Sensors* 21, 4525 (2021).

Jung-Hun Lee, Ji-Soo Kim, Hyun-Joong Kim, Kyujong Park, Jungwoo Moon, Jinyoung Lee and Youngju Park "Free Volume Effect via Various Chemical Structured Monomers on Adhesion Property and Relative Permittivity in Acrylic Pressure Sensitive Adhesives." *Polymers* 12, 2633 (2020) Aleksei Almaev, Vladimir Nikolaev, Pavel Butenko, Sergey Stepanov, Aleksei Pechnikov, Nikita Yakovlev, Igor Sinyugin, Sevastian Shapenkov, Mikhail Scheglov "Gas Sensors Based on Pseudohexagonal Phase of Gallium Oxide." *Phys Status Solidi B Basic Res* 259, 2100306 (2022).

Shuvra Mondal, Bok Ki Min, Yoonsik Yi, and Choon-Gi Choi "Highly Sensitive and Fast Responsive Humidity Sensor based in 2D PtSe2 with Gamma Radiation Tolerance." *Adv. Mater. Technol.* 7, 2100751 (2022).

Mohammad Moein Mohammadi, Abhishek Kumar, Jun Liu, Yang Liu, Thomas Thundat, and Mark T. Swihart "Hydrogen Sensing at Room Temperature Using Flame-Synthesized Palladium-Decorated Crumpled Reduced Graphene Oxide Nanocomposites." *ACS Sens.* 5, 2344-2350 (2020).

Jianwen Sun, Dong Hu, Zewen Liu, Luke M. Middelburg, Sten Vollebregt, Pasqualina M. Sarro, Guoqi Zhang "Low power AlGaN/GaN MEMS pressure sensor for high vacuum application." *Sens. Actuator A Phys.* 314, 112217 (2020).

Seung Bae Son, Ji Won Youn, Kwang-Seok Kim, Dae Up Kim "Optical properties of periodic micropatterned VO_2 thermochromic films prepared by thermal and intense pulsed light sintering." *Mater. Des.* 182, 107970 (2019).

Seyoung Park, Soo-Min Lee, Jin-Kyo Jeong, Donggu Kim, Hyunsoo Kim, Hyun-Sook Lee, Wooyoung Lee "Sensing performance of Pd nanogap supported on an elastomeric substrate in a wide temperature range of -40 to 70 °C." Sens. Actuators B Chem. 348, 130716 (2021).



Osbel Almor, Gebhard J. Matt, Albert These, Andrii Kanak, levgen Levchuk, Shreetu Shrestha, Andres Osvet, Christoph J. Brabec, and Germà Garcia-Belmonte "Surface versus Bulk Currents and Ionic Space-Charge Effects in CsPbBr3 Single Crystals." . *Phys. Chem. Lett.* 13, 824-3830 2022).

Luke M. Middelburg, H. W. van Zeijl, Sten Vollebreg.; Bruno Morana, Guoqi Zhang, "Toward a Self-Sensing Piezoresistive Pressure Sensor for All-SiC Monolithic Integration." *IEEE Sens. J.* 20, 11265-11274 (2020).

Mohd Faizol Abdullah, Nur Julia Nazim Bulya Nazim, Nurhidaya Soriadi, Siti Aishah Mohamad Badaruddin, Mohd Rofei Mat Hussin, Mohd Ismahadi Syono "Thermal Characterization of Multi-Layer Graphene Heat Spreader by Pt/Cu/Ti Micro-Coil." *PHYS STATUS SOLIDI A* 218, 2100301 (2021).

Murat Onen, Nicolas Emond, Ju Li, Bilge Yildiz, and Jesús A. del Alamo "CMOS-Compatible Protonic Programmable Resistor Based on Phosphosilicate Glass Electrolyte for Analog Deep Learning." *Nano Lett.* 21, 6111-6116 (2021).

Shuvra Mondal, Seong Jun Kim, and Choon-Gi Choi "Honeycomb-like MoS₂ Nanotube Array-Based Wearable Sensors for Noninvasive Detection of Human Skin Moisture. " ACS Appl. Mater. Interfaces 12, 17029-17038 (2020).

Jong Won Lee, Jeonghyeon Park, Heera Kwon, Woong-Ki Hong, Jong Kyu Kima and Jaehee Cho "Self-protective GalnN-based light-emitting diodes with VO_2 nanowires." *Nanoscale* 11, 18444-18448 (2019).

Jiarui Mo, Luke M. Middelburg, Bruno Morana, H. W. van Zeijl, Sten Vollebregt, Guoqi Zhang "Surface-Micromachined Silicon Carbide Pirani Gauges for Harsh Environments." *IEEE Sens. J.* 21, 1350-1358 (2020).

Wen Sui, Haoran Wang, Jaesung Lee, Afzaal Qamar, Mina Rais-Zadeh, Philip X.-L. Feng "AlScN-on-SiC Thin Film Micromachined Resonant Transducers Operating in High-Temperature Environment up to 600 °C." Adv. Funct. Mater. 32, 2202204 (2022).

Taejung Kim, Wootaek Cho, Beomsang Kim, Junyeong Yeom, Yeong Min Kwon, Jeong Min Baik, Jae Joon Kim, Heungjoo Shin "Batch Nanofabrication of Suspended Single 1D Nanoheaters for Ultralow-Power Metal Oxide Semiconductor-Based Gas Sensors." *Small* 18, 2204078 (2022).

A.V. Almaev, V.I. Nikolaev, N.N. Yakovlev, P.N. Butenko, S.I. Stepanov, A.I. Pechnikov, M.P. Scheglov, E.V. Chernikov "Hydrogen sensors based on Pt/α -Ga₂O₃:Sn/Pt structures." Sens. Actuators B Chem. 364, 131904 (2022).

Yiyang Li, Elliot J. Fuller, Joshua D. Sugar, Sangmin Yoo, David S. Ashby, Christopher H. Bennett, Robert D. Horton, Michael S. Bartsch, Matthew J. Marinella, Wei D. Lu, A. Alec Talin "Filament-Free Bulk Resistive Memory Enables Deterministic Analogue Switching." Adv. Mater. 32, 2003984 (2020).

Yiyang Li, Elliot J. Fuller, Shiva Asapu, Sapan Agarwal, Tomochika Kurita, J. Joshua Yang, and A. Alec Talin "Low-Voltage, CMOS-Free Synaptic Memory Based on Li_xTiO₂ Redox Transistors." *ACS Appl. Mater. Interfaces* 11, 38982-38992 (2019).





Yiyang Li, Jorik van de Groep, A. Alec Talin, and Mark L. Brongersma "Dynamic Tuning of Gap Plasmon Resonances Using a Solid-State Electrochromic Device." *Nano Lett.* 19, 7988-7995 (2019).

Diana S. Kim, Virgil J. Watkins, Laszlo A. Cline, Jingxian Li, Kai Sun, Joshua D. Sugar, Elliot J. Fuller, A. Alec Talin, Yiyang Li "Nonvolatile Electrochemical Random-Access Memory under Short Circuit." *Adv. Electron. Mater.* 9, 2200958 (2023).

Alexander Stangl, Adeel Riaz, Laetitia Rapenne, José Manuel Caicedo, Juan de Dios Sirvent, Federico Baiutti, Carmen Jiménez, Albert Tarancón, Michel Mermoux and Mónica Burriel "Tailored nano-columnar La₂NiO₄ cathodes for improved electrode performance†‡." *J. Mater. Chem.* A 10,2528-2540 (2022).

Hyeonghun Kim, Woochul Kim, Jiyoon Park, Namsoo Lim, Ryeri Lee, Sung Jun Cho, Yogeenth Kumaresan, Myoung-Kyu Oh and Gun Young Jung, "Surface conversion of ZnO nanorods to ZIF-8 to suppress surface defects for a visible-blind UV photodetectort." *Nanoscale* 10, 21168-21177 (2018).

Y. Pak, Y. Kim, N. Lim, J.-W. Min, W. Park, W. Kim, Y. Jeong, H. Kim, K. Kim, S. Mitra, B. Xin, T.-W. Kim, I. S. Roqan, B. Cho, G.-Y. Jung "Scalable integration of periodically aligned 2D-MoS₂ nanoribbon array." *APL Mater.* 6, 76102 (2018).

Yusin Pak, Yeonggyo Jeong, Naresh Alaal, Hyeonghun Kim, Jeonghoon Chae, Jung-Wook Min, Assa Aravindh Sasikala Devi, Somak Mitra, Da Hoon Lee, Yogeenth Kumaresan, Woojin Park, Tae-Wook Kim, Iman S. Roqan, Gun-Young Jung "Highly Stable and Ultrafast Hydrogen Gas Sensor Based on 15 nm Nanogaps Switching in a Palladium–Gold Nanoribbons Array." *Adv. Mater. Interfaces* 6, 1801442 (2019). M. Rahimi, K. Sobnath, F. Mallet, P. Lafarge, C. Barraud, W. Daney de Marcillac, D. Fournier, and M.L. Della Rocca "Complete Determination of Thermoelectric and Thermal Properties of Supported Few-Layer Two-Dimensional Materials." *Phys. Rev. Appl.* 19, 34075 (2023).

Alexander Stangl, Dolors Pla, Caroline Pirovano, Odette Chaix-Pluchery, Federico Baiutti, Francesco Chiabrera, Albert Tarancón, Carmen Jiménez, Michel Mermoux, Mónica Burriel "Isotope Exchange Raman Spectroscopy (IERS): a novel technique to probe physicochemical processes in situ" Advanced Materials, 2303259 (2023).

Joonchul Shin, Geonhee Lee, Myungwoo Choi, Huiwon Jang, Yunsung Lim, Gwang-Su Kim, Sang-Hyeon Nam, Seung-Hyub Baek, Hyun-Cheol Song, Jihan Kim, Chong-Yun Kang, Jeong-O Lee, Seokwoo Jeon, Donghwi Cho,Ji-Soo Jang "Atomically mixed catalysts on a 3D thin-shell TiO 2 for dual-modal chemical detection and neutralization" *J. Mater. Chem. A.* (2023)

Hyerin Jo, Asad Ali, Won Suk Oh, Sung Jin An, Gyu-Chul Yi, and Hongseok Oh, "Flexible Photonic Synapses Using Vertical ZnO Nanotubes on Graphene Films" *IEEE Sens.* J.30 (2023).

S. Panisset, A. Riaz, A. Stangl, M. Burriel, D. Jauffres, "Contribution of numerical modelling to design oxygen electrode for micro-solid oxide cells: A case study of high-performance nano-columnar La₂NiO₄ thin films" *Journal of Power Sources*, 593: 233951.(2024)

Asad Ali, Jamin Lee, Kyoungho Kim, Hongseok Oh, Gyu-Chul Yi, "Highly Sensitive and Fast Responding Flexible Force Sensors Using ZnO/ZnMgO Coaxial Nanotubes on Graphene Layers for Breath Sensing" *Advanced Healthcare Materials*, 2304140.(2024)



Jong-Hyun Kwak, Wootaek Cho, Beomsang Kim, Taejung Kim, Heungjoo Shin "Hightemperature adaptive through-silicon via with pyrolyzed carbon via-sealing plates for packaging 3D carbon nanostructure-based devices fabricated using C-MEMS" *Carbon, 216,* 118534(2024).

Wootaek Cho, Taejung Kim, Heungjoo Shin "Thermal conductivity detector (TCD)-type gas sensor based on a batch-fabricated 1D nanoheater for ultra-low power consumption" Sensors and Actuators B: Chemical, 371: 132541. (2022)

Shalik Ram Joshi, Beomsang Kim, Shin-Kwan Kim, Wonho Song, Kibog Park, Gun-Ho Kim1 and Heungjoo Shin "Low-Cost and Fast-Response Resistive Humidity Sensor Comprising Biopolymer-Derived Carbon Thin Film and Carbon Microelectrodes" *Journal of Power Journal* of *The Electrochemical Society*, 167.14: 147511.(2020)

Mahesh Peddigari, Bo Wang, Rui Wang, Woon-Ha Yoon, Jongmoon Jang, Hyunjong Lee, Kyung Song, Geon-Tae Hwang, Kai Wang, Yuchen Hou, Haribabu Palneedi, Yongke Yan, Han Seung Choi, Jianjun Wang, Aravindkrishna Talluri, Long-Qing Chen, Shashank Priya, Dae-Yong Jeong, Jungho Ryu "Giant Energy Density via Mechanically Tailored Relaxor Ferroelectric Behavior of PZT Thick Film" *Advanced Materials*, 35.45: 2302554. (2023)

Geon Lee, Tek Kyoung Sung, Hyunseok Song & Jungho Ryu "Electrical properties of $AgNbO_3$ thick films by aerosol deposition followed by post-annealing in oxygen atmosphere" *Journal of the Korean Ceramic Society* 1-10.(2024)

Geon Gug Yang, Dong-Ha Kim, Sanket Samal, Jungwoo Choi, Heejung Roh, Camille E. Cunin, Hyuck Mo Lee, Sang Ouk Kim, Mircea Dincă, Aristide Gumyusenge "Polymer-Based Thermally Stable Chemiresistive Sensor for Real-Time Monitoring of NO₂ Gas Emission"*ACS sensors*, 8.10: 3687-3692.(2023) Arjun Wadhwa, Jaime Benavides-Guerrero, Mathieu Gratuze, Martin Bolduc, and Sylvain G. Cloutier "All Screen Printed and Flexible Silicon Carbide NTC Thermistors for Temperature Sensing Applications" *Materials.* 17(11), 2489 (2024).

Paul Fourmont, Fabrice Vaussenat, Mathieu Gratuze, Caroline A. Ross, Sylvain G. Cloutier "Highly-Sensitive and High Operating Range Fully-Printed Humidity Sensors Based on BiFeO₃/BiOCl Heterojunctions" *Advanced Electronic Materials*, 10(11), 2400156 (2024).

Hyeon-Seo Choi, Minjee Ko, Taejin Lee, Jin-Woo Jung, Young-Jun Lee, Hyeonjong Jeong, Youngjae Kim, Dongha Kim, Jinhee Heo, Shinbuhm Lee, JaeDong Lee, Chang-Hee Cho "Tunable Polariton Rabi Oscillation in Phase-Changing Perovskite Microcavities" Advanced Science, Early view, 2417596 (2025).

Jingxian Li, Andrew J. Jalbert, Sangyong Lee, ..., Elliot J. Fuller, A. Alec Talin, Yiyang Li "Nonvolatile electrochemical memory at 600 °C enabled by composition phase separation" *Device* (2024).





Video Link Sites

1. The Exclusive Smart Probe





2. Nextron Probe's



4. MPS 3D Illustration





5. XYZ Ext. Positioner Instruction





6. MPS Temp. Cycle

Test

3. Probe-sample Test

7. Piezo Module Probe















MEMO





EUROPE



Customers

EUROPE





EUROPE



Customers

AMERICA







ASIA



Customers

KOREA



ETRI

KRİCT

한국화학연구원

KERI

한국전기연구원

圈

서울대학교







dor	acu	14		ø
uUI	991	ary .	2.0	P
LININ	/FRS	ITV.		٠

(laj

Ser all





٠

KŤT

만전생평가면구소.

NINT

나노융합기술원

🕞 LG화학

C LG Electronics

Ø



CANAL U

KITECH 한국생산기술연구표

==(

Safety Solution Lock 쏠락

서강대학교 SOGANG UNIVERSITY







40.4

*NEXTRON

MICRO PROBE SYSTEM







KOREA







NEXTRON

Headquarter

KOREA

609 V1 Tower, 273-20 Gaejwa-Ro, Geumjeong-Gu, Busan-City, 46257, Korea Phone: +82 51 512 6770 Fax: +82 51 512 6737

Marketing contact: sales@nextron.co.kr www.microprobesystem.com

Overseas Branch Office

EUROPE

Paris, France Marketing contact: ykkim@nextron.co.kr **AMERICA**

Chicago, USA Marketing contact: dhkim@nextron.co.kr

Partner

JAPAN

MTT Instruments, Inc., Osaka Marketing contact: t.fukushima@mtt-ins.co.jp

CHINA

www.miniprobesystem.com Marketing contact: 15191589801@163.com